



Science For A Better Life

Jaydess:

Differentiation from Mirena



Aim of this presentation

- Until now only one intrauterine system (IUS) has been approved, Mirena
- Jaydess, a new IUS, is now also approved
- There are differences between these 2 products
- It is important that, in particular, HCPs are informed of how to differentiate between the appearances of Jaydess and Mirena



Duration and indication of use

- Both intrauterine systems, IUS, are T shaped devices that are placed within the uterine cavity and provide a continuous release of levonorgestrel,
- Jaydess is licensed in the indication of contraception for up to 3 years
- Mirena is licensed in the indication of contraception, heavy menstrual bleeding and endometrial protection during use of estrogen replacement therapy, for up to 5 years.



Mirena branding





Jaydess Branding





Visual appearance

28mm horizontal width →



30mm
Vertical
length

Silver
ring

Jaydess

32mm horizontal width →



32mm
vertical
length


Mirena




Jaydess: Ultrasound Images



3D
image







3D
Image



Mirena: Ultrasound images







Science For A Better Life

Thank you!